



Postfach einrichten (Thunderbird)

Sollten Sie Ihr Exchange-Postfach in Thunderbird einbinden wollen, gehen Sie bitte wie folgt vor.


1. Starten Sie Thunderbird und rufen Sie über den Menüpunkt „Datei“ > „Neu“ den Punkt „Bestehendes E-Mail-Konto...“ auf.



2. Im sich öffnenden Fenster geben Sie bitte **Ihren Namen, Ihre E-Mailadresse und Ihr TUB-Kennwort** ein (1). Klicken Sie anschließend auf die Schaltfläche „Manuell einrichten...“ (2).

A screenshot of the 'Bestehende E-Mail-Adresse einrichten' (Existing Email Address Setup) dialog box. The dialog has a title bar with a close button. The main text says 'Bestehende E-Mail-Adresse einrichten' and 'Richten Sie Ihre derzeitige E-Mail-Adresse ein.' Below this are three input fields: 'Ihr Name:' with a placeholder '<Ihr Name>', 'E-Mail-Adresse:' with a placeholder '<Ihre E-Mailadresse>', and 'Passwort:' with a masked password '.....'. There is a checkbox for 'Passwort speichern' (Save password) which is checked. At the bottom, there are three buttons: 'Manuell einrichten...' (highlighted with a red box and a red circle with the number 2), 'Abbrechen' (Cancel), and 'Weiter' (Next). A red box and a red circle with the number 1 highlight the three input fields.

3. Ergänzen Sie im nächsten Schritt die **notwendigen Daten** (1) und klicken Sie anschließend auf „Fertig“ (2), um die Postfacheinrichtung abzuschließen.

A screenshot of the 'Erweiterte Einstellungen' (Advanced Settings) dialog box for an email account. The dialog is divided into two columns: 'POSTEINGANGS-SERVER:' (Incoming Server) and 'POSTAUSGANGS-SERVER:' (Outgoing Server). The 'Protokoll:' (Protocol) is set to IMAP for incoming and SMTP for outgoing. The 'Server:' field contains 'mail.tu-berlin.de' for both. The 'Port:' is 993 for incoming and 587 for outgoing. The 'SSL:' is set to SSL/TLS for incoming and STARTTLS for outgoing. The 'Authentifizierung:' (Authentication) is set to 'Passwort, normal' (Password, normal) for both. The 'Benutzername:' (Username) is '<Kontoname>@tu-berlin.de' for both. At the bottom, there are three buttons: 'Abbrechen' (Cancel), 'Erneut testen' (Test Again), and 'Fertig' (Done, highlighted with a red box and a red circle with the number 2). A red box and a red circle with the number 1 highlight the server configuration fields.

4. Für die Einrichtung benötigte Daten:

	Posteingangs-Server	Postausgangs-Server
Protokoll	IMAP	SMTP
Server	mail.tu-berlin.de	
Port	993	587
SSL	SSL	STARTTLS
Authentifizierung	Passwort, normal	
Benutzername	<Ihr TUB-Accountname>@tu-berlin.de	